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Thành phần:

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Product: RGML PRP Kit

Components:

Spinal anesthesia needle; 5 mL Luer-lock Syringe; Filter; Connector; 10 mL Luer-lock Syringe; Alcohol pad; Kim Vacuum blood collection needle; Needle holder; Vacuum blood collection tube; Bandage; PRP separation tube.

Reference number: 148

Amount: kit/box.

Purpose of use:

For research or manufacturing use only.

RGML PRP Kit kit is used to isolate and activate platelet-rich plasma from 12 of peripheral blood, umbilical cord blood, etc.

Summary of salient features

RGML PRP Kit is a kit for blood collection and manufacturing of activated platelet-rich plasma. This kit draws blood using vacuum tubes and separates platelets by isopyclic centrifugation. PRP is activated through a filtration membrane, resulting in a final activated PRP product that does not contain any blood cells and has a high purity level. RGML PRP Kit was researched and developed by Stem Cell Institute, HCMUS, VNUHCM, Viet Nam. RGML PRP Kit guarantees to be sterile (0 CFU) with endotoxin level ≤ 1 EU/mL.

Reviewed/tested Applications:

RGML PRP Kit provides a fast and convenient method to collect blood and manufacture platelet-rich plasma from peripheral blood or umbilical cord blood, etc.

Storage conditions and shelf life:

Storage and transportation temperature: Room temperature.

Expiry date: following the nearest expiry date of the components.

Procedure:

5. Make sure the product is intact and within the expiry date.
6. **Blood collection**

- 6.1. Attach Vacuum blood collection needle (butterfly type) to the Needle holder.
- 6.2. Clean the puncture site with an alcohol pad, then perform venipuncture.
- 6.3. Advance the Vacuum blood collection tube into the Needle holder. Blood will automatically flow into the tubes.
- 6.4. Collect blood until the tube is full (tube contains 10-12 mL of blood). Gently invert tubes to mix blood with an additive.

7. PRP Production:

- 7.1. Open the Vacuum blood collection tube and PRP separation tube. Then transfer the blood gradually from the Vacuum blood collection tube into the PRP separation tube.
- 7.2. Place the PRP separation tubes containing blood into the centrifuge (place 1 more counterweight tube in a symmetrical position).
- 7.3. Centrifuge the tubes at **1500g for 10 minutes**.
- 7.4. Gently take out the PRP separation tubes from the centrifuge.
- 7.5. Attach the 10 mL Luer-lock Syringe to a Spinal anesthesia needle. Remove supernatant on the top and leave **1,0-1,5 mL of plasma** at the bottom (containing platelet pellet).
- 7.6. Suspend the platelet pellet in the plasma, and aspirate 3-5 times. This mixture is inactivated platelet-rich plasma (PRP).

8. Activation and sterilization of PRP:

- 8.1. Collect inactivated PRP into the 10 mL Luer-lock Syringe.
- 8.2. Attach the 10 mL Luer-lock Syringe to the Filter.
- 8.3. Attach the Connector to the other end of the Filter.

Regenmedlab

Cellatist

Stem food

- 8.4. Attach 5 mL Luer-lock Syringe to the other end of the Connector.
- 8.5. Strongly push the piston of the 10 mL Luer-lock Syringe so that inactivated PRP goes through the Filter and into the 5 mL Luer-lock Syringe (The product in 5 mL Luer-lock Syringe is activated PRP).

Precautions:

Do not use the kit if the packaging is damaged or cracked, or if one of the kit components is abnormal

Troubleshooting:

Not applied

For buying other products, please visit:

<http://biomedmart.org>

For further information, please contact:

contact@sci.edu.vn; sales@sci.edu.vn

Explanation of symbols and warnings

The symbols present on the kit are explained below

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