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Thành phần:

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Product: RGML PRP Kit V1

Components:

Spinal anesthesia needle; 5 mL Luer-lock Syringe; Filter; Connector; 10 mL Luer-lock Syringe; Alcohol pad; Kim Vacuum blood collection needle; Needle holder; Vacuum blood collection tube; Bandage; PRP separation tube.

Reference number: 166

Amount: kit/box.

Purpose of use:

For research or manufacturing use only.

RGML PRP Kit V1 kit is used to isolate and activate platelet-rich plasma from 12 of peripheral blood, umbilical cord blood, etc.

Summary of salient features

RGML PRP Kit V1 is a kit for blood collection and manufacturing of activated platelet-rich plasma. This kit draws blood using vacuum tubes and separates platelets by isopyclic centrifugation. PRP is activated through a filtration membrane, resulting in a final activated PRP product that does not contain any blood cells and has a high purity level. RGML PRP Kit V1 was researched and developed by Stem Cell Institute, HCMUS, VNUHCM, Viet Nam. RGML PRP Kit V1 guarantees to be sterile (0 CFU) with endotoxin level ≤ 1 EU/mL.

Reviewed/tested Applications:

RGML PRP Kit V1 provides a fast and convenient method to collect blood and manufacture platelet-rich plasma from peripheral blood or umbilical cord blood, etc.

Storage conditions and shelf life:

Storage and transportation temperature: Room temperature.

Expiry date: following the nearest expiry date of the components.

Procedure:

5. Make sure the product is intact and within the expiry date.
6. **Blood collection**

- 6.1. Attach Vacuum blood collection needle (butterfly type) to the Needle holder.
- 6.2. Clean the puncture site with an alcohol pad, then perform venipuncture.
- 6.3. Advance the Vacuum blood collection tube into the Needle holder. Blood will automatically flow into the tubes.
- 6.4. Collect blood until the tube is full (tube contains 10-12 mL of blood). Gently invert tubes to mix blood with an additive.

7. PRP Production:

- 7.1. Open the Vacuum blood collection tube and PRP separation tube. Then transfer the blood gradually from the Vacuum blood collection tube into the PRP separation tube.
- 7.2. Place the PRP separation tubes containing blood into the centrifuge (place 1 more counterweight tube in a symmetrical position).
- 7.3. Centrifuge the tubes at **1500g for 10 minutes**.
- 7.4. Gently take out the PRP separation tubes from the centrifuge.
- 7.5. Attach the 10 mL Luer-lock Syringe to a Spinal anesthesia needle. Remove supernatant on the top and leave **1,0-1,5 mL of plasma** at the bottom (containing platelet pellet).
- 7.6. Suspend the platelet pellet in the plasma, and aspirate 3-5 times. This mixture is inactivated platelet-rich plasma (PRP).

8. Activation and sterilization of PRP:

- 8.1. Collect inactivated PRP into the 10 mL Luer-lock Syringe.
- 8.2. Attach the 10 mL Luer-lock Syringe to the Filter.
- 8.3. Attach the Connector to the other end of the Filter.

Regenmedlab

Cellatist

Stem food
From Science to Food

Dermaloka

- 8.4. Attach 5 mL Luer-lock Syringe to the other end of the Connector.
- 8.5. Strongly push the piston of the 10 mL Luer-lock Syringe so that inactivated PRP goes through the Filter and into the 5 mL Luer-lock Syringe (The product in 5 mL Luer-lock Syringe is activated PRP).

Precautions:

Do not use the kit if the packaging is damaged or cracked, or if one of the kit components is abnormal

Troubleshooting:

Not applied

For buying other products, please visit:

<http://biomedmart.org>

For further information, please contact:

contact@sci.edu.vn; sales@sci.edu.vn

Explanation of symbols and warnings

The symbols present on the kit are explained below

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