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SẢN PHẨM/Product: TISSUECRYOSAVE

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*For research or manufacturing use only.*

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*TissueCryosave is a preservation medium for tissues at -86°C (for 12 months) or -196°C (long-term).*

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*TissueCryosave is a complete cryopreservation medium for tissues, requiring no additional components to be added before use.*

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*The formula has been improved by reducing DMSO (10%) and eliminating animal-derived components (xeno-free). In TissueCryosave, bovine serum is replaced with human serum albumin, minimizing the risk of immune responses from cross-contamination of animal-derived components to humans.*

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*TissueCryosave has a low endotoxin level ( $\leq 0.5$  EU/mL), is free from Mycoplasma contamination, meets the sterility requirements (0 CFU), and has a pH range of 6.5-7.4.*

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*TissueCryosave has been evaluated for effective cryopreservation of skin, umbilical cord tissue, and placental tissue. Tissue samples maintained viability after 12 months at 86°C and long-term at -196°C.*

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*Note: The viability of the tissue depends on the state of the tissue at the time of preservation.*

### Hoàn nguyên, pha loãng, phối trộn/ Reconstitution, Dilution, and Mixing:

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*The product is supplied in a 1X form and does not require any dilution or additional components before use.*

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*Transport and storage: -20°C.*

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*Expiry date: 24 months from the date of manufacture..*

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*The vial of TissueCryosave should be thawed overnight prior to use. Thoroughly mix the product after thawing. Ensure that TissueCryosave has reached a temperature of 2-8°C.*

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*The tissue should be cut into small pieces measuring 1-3 mm square and washed with Washing Buffer solution.*

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*To ensure high tissue viability after thawing, tissue fragments should be sharply cut, blood should be thoroughly washed off, and necrotic areas should be removed.*

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*Gently transfer 20-50 tissue fragments into a 2 mL cryovial*

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*Handle with care, avoiding tissue damage.*

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*Draw TissueCryosave solution into the cryovial, filling approximately 80% of the cryovial's capacity. Securely tighten the cryovial cap.*

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*Cool the cell containers by placing them in Mr. Frosty or using a programmed cooling device set at a rate of -1°C per minute until reaching -86°C. After reaching -86°C, store the samples at -86°C overnight (before transferring to -196°C if storing at -196°C)*

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*To ensure a high cell viability after thawing, the time from contact with TissueCryosave to cooling should not exceed 15 minutes*

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*The cooling process using Mr. Frosty is carried out by placing the cell-containing tubes into the box and placing it in a minimum of -86°C freezer for at least 4 hours or overnight before transferring the cell tubes to another appropriate storage container. Afterward, the Mr. Frosty should be removed from the freezer and kept at room temperature when not in use (before and after cooling). It is recommended to replace the isopropanol every 5 freezing cycles*

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*Cells can be stored at either -86°C or -196°C, depending on the application and the specific production process.*

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Remove the tissue-containing tubes/bags from the -86°C freezer or liquid nitrogen and place them in a 37°C incubation bath.

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Wait for the tubes/bags to thaw completely. If there is still a small amount of ice inside the tubes, promptly transfer them to a biosafety cabinet

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Note: DMSO in the TissueCryosave becomes highly toxic to cells at temperatures above 4 degrees Celsius; therefore, after the ice has just melted, quickly rinse and remove the cryopreservation solution.

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Wash the tissue twice with Washing Buffer.

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The results will be better if Thawbest is used to wash the tissue after thawing.

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Note: Cell growth from tissue may be slower compared to fresh tissue culture. In case cells do not grow from the tissue and tissue fragments do not adhere to the culture surface, proceed by cutting the tissue fragments into smaller pieces and re-culturing under the same conditions.

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TissueCryosave is not a medicine and is not used as a drug or infusion solution.

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Do not use this product if the packaging is damaged or cracked, and/or if the solution appears discolored or turbid.

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The symbols present on the kit are explained below:

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